


Achalasia Balloon Dilator

Intended for dilation of the lower esophageal sphincter in patient with achalasia.

Technical Information:

- Inflated Balloon Diameter - 25mm, 30mm, 35 mm
- Balloon length : 80 mm.
- Shaft Diameter : 5mm.
- Working Length : 90 cm.
- Wire Guide Compatibility : 0.038"
- Burst pressure : 3 ATM
- High Radial Vector Force.
- Atraumatic Tip.
- Radioopaque marker inside the balloon.
- Sterile for Single use only.



Ordering Information-Mention product code x quantity required

Reference	Product Name	BalloonØ
SABD - 25	Achalasia Balloon Dilator	25 mm
SABD - 30	Achalasia Balloon Dilator	30 mm
SABD - 35	Achalasia Balloon Dilator	35 mm