

Biopsy Forcep (Gastro)

Used endoscopically to obtain mucosal tissue samples from bronchiole for microscopic histopathological examination.

Technical Information:

- Availability - Sterile For single use only.
- Compatible With - 2.0 / 2.8 mm instrument channel
- Length : 110 Cm
- Individually Packed



Ordering Information-Mention product code x quantity required

Reference	Product Name	Type	Required ChannelØ	Working Length
S- BF- WS-17	Biopsy Forcep	With Spike	2.8 mm	170cm
S-BF-WS-21	Biopsy Forcep	Without Spike	2.8 mm	210cm
S- BF-17	Biopsy Forcep	With Spike	2.8 mm	170cm
S-BF-21	Biopsy Forcep	Without Spike	2.8 mm	210cm