


Cleaning Brush (Gastro)

Used to clean endoscope channel.

Technical Information:

- Type both end bristles
- Shaft plastic coated wire
- Semi soft and channel friendly bristles for cleaning of endoscope channel.
- Atraumatic tip
- Compatible with 2.0 / 2.8 mm instrument channel
- Length 110cm
- Individually packed



Ordering Information-Mention product code x quantity required

Reference	Product Name	Required Channel Ø	Working Length
S - CB - DE	Cleaning Brush Double Ended	2.8 mm.	230 cm
S - CB - DE-BS	Cleaning Brush Double Ended Big-Small	2.8 mm.	230 cm
S - CB - P	Port Cleaning Brush	-	-
S - CB - B	Port Cleaning Brush	-	-
S - CB - RG	Cleaning Brush Reusable Gastro	2.8 mm	180cm
S - CB - RC	Cleaning Brush Reusable Colono	2.8 mm	230cm
S - CB - D	Cleaning Brush for Dilator	1.2 mm	120cm