

Clips for Endoclip


Used to mechanically bind blood vessels

Indications

- Hemostasis for: mucosal/submucosal
- defects <3cm, bleeding ulcers/arteries <2mm,
- surgery sites, closure of GI luminal performance

Technical Information

- Jaw angle=1350
- Distance between open clips>8mm
- Clips for endoclip is supplied sterile



Ordering Information-Mention product code x quantity required

S-CL-EC Minimum order Qty. is 20pcs/box

Reference	Product Name	Width W	Distance Between open Clips H	JAW angle A
S-CL-EC	Clips for Endoclips	1.4mm	≥ 8mm	1350