


Endoknot

Used for Endoscopic knotting before cut / coagulation in GI to prevent excessive bleeding.

Bio absorbable thread enhance more indications

Technical Information

- Working length 230cm
- Compatible with instrument channel 2.8mm
- Loop Diameter 30mm & 40mm
- Sheath diameter 2.5 mm
- Supplied sterile and it is for single use



Ordering Information-Mention product code x quantity required

Reference	Product Name	Loop Ø	Required ChannelØ	Working Length
S-EK-30	ENDOKNOT	30mm	2.8mm	230cm
S-EK-40	ENDOKNOT	40mm	2.8mm	230cm