


Esophageal Gastric Balloon Retrieval Device

Used for retrieval of Esophageal Gastric Balloon

Technical Information

- Compatible with 2.8 mm. instrument channel.
- Sheath diameter 2.5 mm
- Working length 170 cm.
- Supplied Sterile, for single use only .
- Excellent and secure grip to retrieve deflated Esophageal Gastric Balloon



Ordering Information-Mention product code x quantity required

Reference	Product Name	Working Length	Required ChannelØ
SEGBRD	Esophageal Gastric Balloon Retrieval Device	170 cm	2.8 mm