


Foreign Body Holding Forcep (2,3,4 Prong) (Gastro)

Used to retrieve Foreign body from G.I. Tract.

Technical Information:

- Compatible with 2.0 & 2.8 mm. instrument channel.
- Working length 170 cm. / 230 cm.
- Supplied Non Sterile and It is Reusable.
- Excellent and secure grip to retrieve Foreign
- Body from G.I. Tract.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Working Length	Required Channel Ø	Jaws
SFBH-223	F B HOLDER 2 PRONG	170 cm	2.8 mm	2 Prong
SFBH-218	F B HOLDER 2 PRONG	170 cm	2.0 mm	2 Prong
SFBH-323	F B HOLDER 3 PRONG	170 cm	2.8 mm	3 Prong
SFBH-318	F B HOLDER 3 PRONG	170 cm	2.0 mm	3 Prong
SFBH-423	F B HOLDER 4 PRONG	170 cm	2.8 mm	4 Prong
SFBH-223-C	F B HOLDER 2 PRONG	230 cm	2.8 mm	2 Prong
SFBH-218-C	F B HOLDER 2 PRONG	230 cm	2.0 mm	2 Prong
SFBH-323-C	F B HOLDER 3 PRONG	230 cm	2.8 mm	3 Prong
SFBH-423-C	F B HOLDER 4 PRONG	230 cm	2.8 mm	4 Prong

“Shaili “ Range of Foreign Body Holder : 2 prong , 3 Prong , Retrieval Basket, Snare, Shark Tooth Forcep, Alligator Forcep, Monofilament Memory Basket, Loop Basket.