


Hot Biopsy Forcep

Used endoscopically to obtain mucosal tissue samples from Gastrointestinal tract for microscopic- Histopathological Examination.

Monopolar electro surgical unit can be used to cauterize the tissue when bleeding occurs during tissue sampling

Technical Information:

- Without Spike
- Availability- Sterile, for single use
- Compatible with 2.8 instrument channel
- Working Length- 170, 210cm
- Active cord available separately.
- Individually packed



Ordering Information-Mention product code x quantity required

| Reference | Product Name | Type | Required ChannelØ | Working Length |
|-----------|-------------------|---------------|-------------------|----------------|
| S-BF-H-17 | Hot Biopsy Forcep | Without Spike | 2.8 mm | 170cm |
| S-BF-H-21 | Hot Biopsy Forcep | Without Spike | 2.8 mm | 210cm |