


Litho Basket-TTS

Used for purpose of CBD stone retrieval/ lithotripsy of stone in bile duct through the scope

Technical Information:

- Working length 230cm
- Extension wire length 205cm
- Loop Diameter 25/30/35/40/45/50/55mm
- Sheath Diameter 2.5mm
- Recommended for manual lithotripsy.
- Supplied sterile and it is for single use
- Detachable Handle.



Ordering Information-Mention the product code x quantity required.

Reference	Product Name	Basket Ø	SheathØ	Working Length
SLTB-25	Lithotripsy Basket TTS	25mm	2.5mm	230cm
SLTB-30	Lithotripsy Basket TTS	30mm	2.5mm	230cm
SLTB-35	Lithotripsy Basket TTS	35mm	2.5mm	230cm
SLTB-40	Lithotripsy Basket TTS	40mm	2.5mm	230cm
SLTB-45	Lithotripsy Basket TTS	45mm	2.5mm	230cm
SLTB-50	Lithotripsy Basket TTS	50mm	2.5mm	230cm
SLTB-55	Lithotripsy Basket TTS	55mm	2.5mm	230cm