


Loop Basket

Used for Endoscopic retrieval of Foreign Bodies from the GI Tract and to retrieve Polyp through Colono Scope.

Technical Information:

- Shape : Oval
- With flush port.
- Availability - Non sterile and for single use only.
- Loop Diameter > 61 mm could be with 2 hand operation handle



Ordering Information-Mention product code x quantity required

Reference	Product Name	Required ChannelØ	Loop Size	Length
SFBH-L-30	LOOP BASKET	2.8mm	30mm	170cm
SFBH-L-40	LOOP BASKET	2.8mm	40mm	170cm
SFBH-L-60	LOOP BASKET	2.8mm	60mm	170cm
SFBH-L-90	LOOP BASKET	2.8mm	90mm	170cm
SFBH-L-C-30	LOOP BASKET	2.8mm	30mm	230cm
SFBH-L-C-40	LOOP BASKET	2.8mm	40mm	230cm
SFBH-L-C-60	LOOP BASKET	2.8mm	60mm	230cm
SFBH-L-C-90	LOOP BASKET	2.8mm	90mm	230cm
SFBH-LB-DBE-60	LOOP BASKET	2.8mm	60mm	280cm