


Memory Basket

Used to retrieve foreign body / Stone from G.I.Tract - through the scope.

Technical Information:

- Seven basket sizes to facilitate foreign body retrieval from G.I.Tract / stone retrieval from bileduct.
- Basket Type - Oval, Hexagonal, Tipless & Helical
- Retains shape for longer time.
- Handle - Light weight with flush port for contrast.
- **Compatible with Lithotripter.**
- Availability - Sterile - for single use only



Ordering Information-Mention product code x quantity required

Reference	Product Name	Required Channel Ø	Basket Ø	Working Length
SMMB-15	Memory Basket	2.8 mm.	15mm	230cm
SMMB-20	Memory Basket	2.8 mm.	20mm	230cm
SMMB-25	Memory Basket	2.8 mm.	25mm	230cm
SMMB-30	Memory Basket	2.8 mm.	30mm	230cm
SMMB-35	Memory Basket	2.8 mm.	35mm	230cm
SMMB-40	Memory Basket	2.8 mm.	40mm	230cm
SMMB-45	Memory Basket	2.8 mm.	45mm	230cm
SMMB-15-T	Memory Tipless Basket	2.8 mm.	15mm	230cm
SMMB-20-T	Memory Tipless Basket	2.8 mm.	20mm	230cm
SMMB-25-T	Memory Tipless Basket	2.8 mm.	25mm	230cm
SMMB-30-T	Memory Tipless Basket	2.8 mm.	30mm	230cm
SMMB-35-T	Memory Tipless Basket	2.8 mm.	35mm	230cm
SMMB-40-T	Memory Tipless Basket	2.8 mm.	40mm	230cm
SMMB-45-T	Memory Tipless Basket	2.8 mm.	45mm	230cm

Tailormade device can be produced provided by pre approved drawing/ specification