

Multi Band Ligator Colonic

Used for Endoscopic Ligation of Colonic Varices.
This device is supplied Non sterile,individually packed-intended for single use only.

- **Technical Information:** Shaili Band Ligator Colonic requires 12.8 mm to 13.5 mm endoscope diameter.
- **Band Ligator Colonic Set Contains** Applicator x 1 pc., Barrel loaded with 6 x Bands,Barrel loaded with 7 x Bands,Atraumatic plastic hook with Tube X1pc.,Flushing Needle blunt tip x1 pc. Second last band with different colour.
- **Barrel Colonic Contains** Barrel loaded with 6 x Bands,Barrel loaded with 7 x Bands,Atraumatic plastic hook with Tube X1pc.



Ordering Information-Mention product code x quantity required

| Reference | Product Name |
|-----------|-----------------------------------|
| SBLS-6-C | MULTI BAND LIGATOR COLONIC SET |
| SBL-6-C | MULTI BAND LIGATOR COLONIC BARREL |
| SBLS-7-C | MULTI BAND LIGATOR COLONIC SET |
| SBL-7-C | MULTI BAND LIGATOR COLONIC BARREL |