

Multi Band Ligator

Used for Endoscopic Ligation of Esophageal Varices at or above Gastro Esophageal Junction. This device is supplied Non sterile, individually packed - intended for single use only.

Technical Information:

- Shaili Band Ligator requires 9.5 to 11.2 .mm endoscope diameter.
- Shaili Paediatric Band Ligator requires 8.0 to 9.5 mm endoscope diameter.

Band Ligator Colonic Set Contains

- Applicator x 1 pc., Barrel loaded with 6 x Bands,
- Barrel loaded with 7 x Bands,
- Atraumatic plastic hook with Tube X1pc.,
- Flushing Needle blunt tip x1 pc.
- Second last band with different colour.

Barrel Colonic Contains

- Barrel loaded with 6 x Bands,
- Barrel loaded with 7 x Bands,
- Atraumatic plastic hook with Tube X1pc.



Ordering Information- Mention the product code x quantity required.

Reference	Product Name	Reference	Product Name
SBLS-6-F*	MULTI BAND LIGATOR SET	SBLS-6-O*	MULTI BAND LIGATOR SET
SBLS-7-F*	MULTI BAND LIGATOR SET	SBLS-7-O*	MULTI BAND LIGATOR SET
SBL-6-F*	MULTI BAND LIGATOR BARREL	SBL-6-O*	MULTI BAND LIGATOR BARREL
SBL-7-F*	MULTI BAND LIGATOR BARREL	SBL-7-O*	MULTI BAND LIGATOR BARREL
SBLS-6-P-F*	MULTI BAND LIGATOR PAEDIATRIC	SBLS-6-P-O*	MULTI BAND LIGATOR PAEDIATRIC
SBL-6-P-F*	MULTI BAND LIGATOR BARREL PAEDIATRIC	SBL-6-P-O*	MULTI BAND LIGATOR BARREL PAEDIATRIC

*To Indicate the Length of Thread, F = 175 cm. for Fujinon Endoscope, O = 155 cm. For Olympus / Pentax Endoscope.