

Peg Kit

Used for Percutaneous Endoscopic Gastrostomy placement To provide enteral nutrition to the patient requiring nutritional support.

Technical Information: Size - 20 -22- 24 Fr. Pull technique Other size can be produced on request. Availability - Sterile-for single use only.

PEG KIT CONTAINS:

- 1 Silicon feeding tube.
- 1 Loop insertion wire,1 Cable tie,1 Universal adaptor,
- 1 Bolus adaptor,1Silicon bush for locking,
- Introducer needle - Big needle,11 No. scalpel,5 ml. Syringe with needle,
- 25 Gauge needle - Small needle,1 Scissor,1 Instruction manual.
- Availability of snare is optional.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Size	Remarks
S-PEG-20	PEG KIT	20FR	
S-PEG-20-SN	PEG KIT	20FR	With Snare
S-PEG-22	PEG KIT	22FR	
S-PEG-22-SN	PEG KIT	22FR	With Snare
S-PEG-24	PEG KIT	24FR	
S-PEG-24-SN	PEG KIT	24FR	With Snare