


Polypectomy Snare

Used for Endoscopic Excision of polyps using monopolar electro-surgical unit and retrieval of foreign bodies from G.I. tract.

Technical Information:

- Availability - Sterile for Single use only.
- With flush port to clean wire during procedure.
- Outer sheath - PTFE.
- Loop Diameter - 6 - 10 - 13 - 15 - 20 - 25 - 30 - 35 - 40 - 45 - 50 - 55 mm
- Loop Diameter for Duckbill Snare is 25 mm only
- Shape - Hexagonal / Oval / Crescent / Duckbill available in Multifilament
- Monofilament available in only Oval shape
- Individually packed .



Ordering Information-Mention product code x quantity required

Reference	Product Name	Length	Shape	Required ChannelØ
SPSO - Loop Ø	Polypectomy Snare	220 cm	OVAL	2.8 mm
SPSH - Loop Ø	Polypectomy Snare	220 cm	HEXAGONAL	2.8 mm
SPSC - Loop Ø	Polypectomy Snare	220 cm	CRESCENT	2.8 mm
SPSD - 25	Polypectomy Snare	220 cm	DUCKBILL	2.8 mm
SPSM - Loop Ø	Polypectomy Snare	220 cm	MONOFIAMENT	2.8 mm
SPSO - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	OVAL	2.8 mm
SPSO - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	HEXAGONAL	2.8 mm
SPSO - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	CRESCENT	2.8 mm
SPSO - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	MONOFILAMENT	2.8 mm
SPSO - Loop Ø-R	Polypectomy Snare Rotatable	220 cm	OVAL	2.8 mm

Tailormade device can be produced provided by pre approved drawing/ specification