

Retrieval Basket


Used to retrieve Foreign Bodies and C.B.D. stones endoscopically

Technical Information:

- Compatible with instrument channel 2.8 mm.
- Basket diameter - 15 mm. / 20 mm. / 25 mm. / 30 mm. / 35 mm. / 40 mm. / 45 mm. / 50 mm.
- Multifilament, Hexagonal shaped.
- Flush port for contrast.
- Atraumatic basket tip for ease of insertion and procedure.

Compatible with Lithotripter.

- Sterile for single use only.



Ordering Information-Strand diameter is 0.41mm, Working length 230 cm.

Reference	Product Name	Basket Ø	Required Channel Ø	Working Length
SRB-4152	Retrieval Basket	15 mm.	2.8 mm.	230 cm.
SRB-4202	Retrieval Basket	20 mm.	2.8 mm.	230 cm.
SRB-4252	Retrieval Basket	25 mm.	2.8 mm.	230 cm.
SRB-4302	Retrieval Basket	30 mm.	2.8 mm.	230 cm.
SRB-4352	Retrieval Basket	35 mm.	2.8 mm.	230 cm.
SRB-4402	Retrieval Basket	40 mm.	2.8 mm.	230 cm.
SRB-4452	Retrieval Basket	45 mm.	2.8 mm.	230 cm.
SRB-4502	Retrieval Basket	50 mm.	2.8 mm.	230 cm.
SRB-4302-DBE	Retrieval Basket	30 mm.	2.8 mm.	280 cm.

Tailormade device can be produced provided by pre approved drawing/ specification