


Sclerotherapy Needle

Used for Endoscopic injection to introduce a sclerosing agent of vasoconstrictor into selected sites to control actual or potential bleeding lesions. The injection of saline to aid in Endoscopic Mucosal Resection (EMR), Polypectomy procedures and to control non-variceal haemorrhage.

Technical specifications:

- Needle - Short and Sharp Bevel, Available in 19/20/21/22/23/24/25 SWG.
- Working length - 170 cm. / 230 cm. / 280 cm
- Availability - Sterile For Single use only.
- Specially Developed Needle To Provide Secure Firm Grip With Inner Tube & Prevent Possible Leakage From The Joint Of Inner Tube & Needle.
- Specially Developed Needle Give "JET" Pressure To Inject The Drug for 23/24/25 SWG.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Needle Gauge	Working Length	Required Channel Ø
SGN-2319	Sclerotherapy Needle	19G	170 cm	2.8 mm.
SGN-2319-C	Sclerotherapy Needle	19G	230 cm	2.8 mm.
SGN-2320	Sclerotherapy Needle	20G	170 cm	2.8 mm.
SGN-2320-C	Sclerotherapy Needle	20G	230 cm	2.8 mm.
SGN-2321	Sclerotherapy Needle	21G	170 cm	2.8 mm.
SGN-2321-C	Sclerotherapy Needle	21G	230 cm	2.8 mm.
SGN-2321-M	Sclerotherapy Needle Metal Tip	21G	170 cm	2.8 mm.
SGN-2321-C-M	Sclerotherapy Needle Metal Tip	21G	230 cm	2.8 mm.
SGN-2322	Sclerotherapy Needle	22G	170 cm	2.8 mm.
SGN-2322-C	Sclerotherapy Needle	22G	230 cm	2.8 mm.
SGN-2323	Sclerotherapy Needle	23G	170 cm	2.8 mm.
SGN-2323-C	Sclerotherapy Needle	23G	230 cm	2.8 mm.
SGN-2323-M	Sclerotherapy Needle Metal Tip	23G	170 cm	2.8 mm.

SGN-2323-C-M	Sclerotherapy Needle Metal Tip	23G	230 cm	2.8 mm.
SGN-2324	Sclerotherapy Needle	24G	170 cm	2.8 mm.
SGN-2324-C	Sclerotherapy Needle	24G	230 cm	2.8 mm.
SGN-2325	Sclerotherapy Needle	25G	170 cm	2.8 mm.
SGN-2325-C	Sclerotherapy Needle	25G	230 cm	2.8 mm.
SGN-2324-DBE	Sclerotherapy Needle	24G	280 cm	2.8 mm.