

Spring Tip Steel Guide Wire

Used to Guide Esophageal Dilators for Dilatation.

- Specially designed Flexible atraumatic tip
- Spring base and metallic wire shaft are interlocked.
- Joint of Spring and wire shaft is tapered to stop guide wire getting stuck in the dilator.
- Smooth surface finish offers easy passage through dilator lumen during dilatation procedure.
- Availability - Sterile, Single Use.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Length	Tip Ø
SMGW	Spring Tip Steel Guide Wire	250 cm	2.00 mm
SMGW-M	Spring Tip Steel Guide Wire Marked	250 cm	2.00 mm