


Stiff guide wire (Gastro)

Used to gain access endoscopically / percutaneous through the difficult strictures.

Technical Information

- Size Diameter - 0.035 inch.
- Availability - Individually packed. Sterile For single use only.
- Length - 110 cm, & 260 cm
- Tip - Straight



Ordering Information-Mention product code x quantity required

Reference	Product Name
S-SGW-035-S-110	Stiff Guide Wire 0.035 Straight
S-SGW-035-S-260	Stiff Guide Wire 0.035 Straight

Tailormade device can be produced provided by pre approved drawing/ specification