


Super Stiff Hydra Neg Guide Wire

Used to gain access through difficult biliary and pancreatic strictures . To facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts.

Technical Information:

- Size Diameter - 0.025 inch./0.032 inch./0.035 inch. / 0.038 inch.
- Length - 450 cm.
- Tip - Straight, J, Angle
- Availability - Individually packed. Sterile For single use only.
- Stiffness of 0.035 with 0.025 size.



Ordering Information-Mention the product code x quantity required.

Reference	Product Name	Wire Ø	Tip
S-SSHYNEG-025-S - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.025 "	STRAIGHT
S-SSHYNEG-032-S - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.032 "	STRAIGHT
S-SSHYNEG-035-S - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.035 "	STRAIGHT
S-SSHYNEG-038-S - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.038 "	STRAIGHT
S-SSHYNEG-025-A - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.025 "	ANGLE
S-SSHYNEG-032-A - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.032 "	ANGLE
S-SSHYNEG-035-A - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.035 "	ANGLE
S-SSHYNEG-038-A - 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.038 "	ANGLE
S-SSHYNEG-025-J- 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.025 "	J Tip
S-SSHYNEG-032-J- 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.032 "	J Tip
S-SSHYNEG-035-J- 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.035 "	J Tip
S-SSHYNEG-038-J- 450	SUPER STIFF HYDRA NEG GUIDE WIRE	0.038 "	J Tip