

## Wire Guided Memory Basket


Used to retrieve foreign body / Stone from G.I.Tract - through the scope.

### Technical Information:

- Three basket sizes to facilitate stone retrieval from biliary tract.
- Basket shapes - Oval, Hexagonal.
- Retains shape for longer time.
- Handle - Light weight with flush port for contrast.
- Centre port for guide wire
- Wire guide compatibility 0.035"

### Not recommended for Lithotripsy.

- Availability - Sterile - for single use only



Ordering Information-Mention product code x quantity required

Reference	Product Name	Required Channel Ø	BasketØ	Working Length
SMMB-WG-30	Memory Basket Wire Guided	2.8mm	30mm	230cm
SMMB-WG-35	Memory Basket Wire Guided	2.8mm	35mm	230cm
SMMB-WG-40	Memory Basket Wire Guided	2.8mm	40mm	230cm