


Endo Lithotripter

Used for Purpose Of C.B.D. Stone Retrieval through the Scope. And for Lithotripsy to Crush Stones in Biliary Duct.

Technical Information:

- Working length- 210 cm.
- Compatible with instrument
- Channel- 4.2 mm.
- Basket shape-Hexagonal
- Flush port for contrast.
- Atraumatic bullet type basket tip for ease of insertion and procedure.



Ordering Information-Mention product code x quantity required

Reference	Product Name	SheathØ	BasketØ	Working Length
SEL-SHEATH	Endolithotripter Flexible Sheath With Basket Filament & Assembly	3.7 mm.	35 mm.	210 cm
SEL - B - 35	Endolithotripter Basket Only		35 mm.	210 cm
SEL - B - 40	Endolithotripter Basket Only		40 mm.	210 cm
SEL - H	Endolithotripter Basket Only			
SLITH-HRB-35	Endolithotripter Basket with PTFE tube		35 mm.	210 cm
SLITH-HRB-40	Endolithotripter Basket with PTFE tube		40 mm.	210 cm