

Esophageal Metal Stent

Used to treat Tracheoesophageal Fistulas and expansion treatment of the Esophageal Stenosis caused by Esophageal Carcinoma, Chemical Injuries, Cardiac Carcinoma, Anastomatic stenosis, and other Traumas, which are surgical contraindicated.

Technical Information

- Radiotherapy compatible.
- Esophageal Stents with Introducer system.
- Expanded Stent diameter. - 18mm.
- Stent wire diameter. - 0.20mm.
- Expanded Stent length - 50mm / 100mm / 120mm / 150mm covered stent.
- Compatible with 0.038 inch wire guide.
- Availability - Sterile - For single use only.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Stent Length
S - ES - 5	Esophageal Metal Stent	50mm
S - ES - 10	Esophageal Metal Stent	100mm
S - ES - 12	Esophageal Metal Stent	120mm
S - ES - 15	Esophageal Metal Stent	150mm